



Anti-PGA polyclonal antibody (CABT-BL2886)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Pepsin from pig stomach.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Pig
Purification	IgG fraction
Conjugate	Unconjugated
Applications	Dot Blot, WB, ELISA, IP, Conjugation
Cellular Localization	Secreted
Format	Liquid
Buffer	0.15M Sodium chloride, 0.02M Potassium phosphate. pH 7.2
Preservative	0.01% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent.
---------------------	--

GENE INFORMATION

Entrez Gene ID [396892](#)

Protein Refseq [NP_999038](#)

UniProt ID [A1YM35](#)
