



Anti-PGCP (aa 71-120) polyclonal antibody (CABT-BL2887)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 71-120 (VDTVGPRLSG SKNLEKAIQI MYQNLQQDGL EKVHLEPVRI PHWERGEESA) of Human PGCP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes a metallopeptidase that belongs to the peptidase M28 family. The encoded protein may catalyze the cleavage of dipeptides with unsubstituted terminals into amino acids.
---------------------	---

GENE INFORMATION

Entrez Gene ID [10404](#)

Protein Refseq [NP_057218](#)

UniProt ID [A0A024R9B8](#)
