



Anti-PHF13 (aa 35-84) polyclonal antibody (CABT-BL2888)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 35-84 (LAYAGYIPYPKEELPLRSSPSPANSTAGTIDSDGWDAGF SDIASSVPLPV) of Human PHF13
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus. Nucleus; nucleoplasm. Predominantly bound to chromatin, but a minor proportion is also detected in the nucleoplasm.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Modulates chromatin structure. Required for normal chromosome condensation during the early
---------------------	---

stages of mitosis. Required for normal chromosome separation during mitosis.

GENE INFORMATION

Entrez Gene ID	148479
----------------	------------------------

Protein Refseq	NP_722519
----------------	---------------------------

UniProt ID	Q86YI8
------------	------------------------
