



# Anti-PHF3 (C-terminal) polyclonal antibody (CABT-BL2890)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide (Human) from C-terminal
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P
<b>Format</b>	Liquid
<b>Buffer</b>	50% Glycerol, 0.88% Sodium chloride, 49.1% PBS pH: 7.40 Note: PBS is without Mg2+ and Ca2+.
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C. Stable for 12 months at -20°C

## BACKGROUND

**Introduction** PHD finger protein 3 is a protein that in humans is encoded by the PHF3 gene.

## GENE INFORMATION

Entrez Gene ID [23469](#)

---

Protein Refseq [NP\\_055968](#)

---

UniProt ID [B3KP41](#)

---