



# Anti-PHLDA3 (aa 38-87) polyclonal antibody (CABT-BL2893)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within N terminal amino acids 38-87 (LQLFEAKGTG GRPKELSFAR IKAVECVEST GRHIYFTLVT EGGGEIDFRC) of Human PHLDA3.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Cytoplasm. Membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	PHLDA3 (pleckstrin homology-like domain, family A, member 3) is a protein-coding gene. Diseases associated with PHLDA3 include germ cell tumors, and lung cancer, and among its
---------------------	---

related super-pathways are Cell Cycle / Checkpoint Control. GO annotations related to this gene include phosphatidylinositol-4,5-bisphosphate binding and phosphatidylinositol-3-phosphate binding. An important paralog of this gene is PHLDA1.

---

## GENE INFORMATION

Entrez Gene ID	<a href="#">23612</a>
----------------	-----------------------

---

Protein Refseq	<a href="#">NP_036528</a>
----------------	---------------------------

---

UniProt ID	<a href="#">Q9Y5J5</a>
------------	------------------------

---