



Anti-PHLDA3 (aa 115-127) polyclonal antibody (CABT-BL2894)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: RARQSLGTGTLVS, corresponding to C terminal amino acids 115-127 of Human PHLDA3.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm. Membrane.
Format	Liquid
Size	100 µg
Buffer	0.5% BSA, Tris buffered saline. pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	PHLDA3 (pleckstrin homology-like domain, family A, member 3) is a protein-coding gene.
---------------------	--

Diseases associated with PHLDA3 include germ cell tumors, and lung cancer, and among its related super-pathways are Cell Cycle / Checkpoint Control. GO annotations related to this gene include phosphatidylinositol-4,5-bisphosphate binding and phosphatidylinositol-3-phosphate binding. An important paralog of this gene is PHLDA1.

GENE INFORMATION

Entrez Gene ID	23612
----------------	-----------------------

Protein Refseq	NP_036528
----------------	---------------------------

UniProt ID	Q9Y5J5
------------	------------------------
