



Anti-PHLDB2 polyclonal antibody (CABT-BL2897)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human PHLDB2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Seems to be involved in the assembly of the postsynaptic apparatus. May play a role in acetylcholine receptor (AChR) aggregation in the postsynaptic membrane.
--------------	--

GENE INFORMATION

Entrez Gene ID [90102](#)

Protein Refseq [NP_665696.1](#)

UniProt ID [Q86SQ0](#)
