



Anti-PIWI (aa 350-450) polyclonal antibody (CABT-BL2905)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from residues 350 - 450 of Drosophila melanogaster Piwi.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Size	100 µg
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [34521](#)

Protein Refseq [NP_476875](#)

UniProt ID [Q9VKM1](#)
