



Anti-PKDREJ polyclonal antibody (CABT-BL2911)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human PKDREJ.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This intronless gene encodes a member of the polycystin protein family. The encoded protein contains 11 transmembrane domains, a receptor for egg jelly (REJ) domain, a G-protein-coupled receptor proteolytic site (GPS) domain, and a polycystin-1, lipoxigenase, alpha-toxin (PLAT) domain. This protein may play a role in human reproduction. Alternative splice variants
--------------	--

have been described but their biological natures have not been determined.

GENE INFORMATION

Entrez Gene ID	10343
----------------	-----------------------

Protein Refseq	NP_006062
----------------	---------------------------

UniProt ID	Q9NTG1
------------	------------------------
