



Anti-PLA2G4E (internal region) polyclonal antibody (CABT-BL2919)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from an internal region of Human PLA2G4E.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Cytoplasm, cytosol. Lysosome membrane, Peripheral membrane protein.
Format	Liquid
Buffer	PBS, 0.88% Sodium chloride, 50% Glycerol pH: 7.40 Note: without Mg ²⁺ and Ca ²⁺
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID	123745
Protein Refseq	NP_001193599

UniProt ID

[B7WPN2](#)
