



Anti-PLXND1 (aa 1800-1900) polyclonal antibody (CABT-BL2930)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1800 - 1900 of Human Plexin D1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane.
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	PLXND1 (plexin D1) is a protein-coding gene. Diseases associated with PLXND1 include ischemic retinopathy, and mobius syndrome. GO annotations related to this gene include
---------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

semaphorin receptor activity. An important paralog of this gene is PLXNA1.

GENE INFORMATION

Entrez Gene ID	23129
----------------	-----------------------

Protein Refseq	NP_055918
----------------	---------------------------

UniProt ID	Q9Y4D7
------------	------------------------
