



Anti-PNMA2 polyclonal antibody (CABT-BL2938)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from the C terminal region of human Paraneoplastic antigen Ma2, conjugated to an immunogenic carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Cellular Localization	Nucleus; nucleolus
Format	Liquid
Buffer	Whole serum, 1X PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C (add glycerol to a final volume of 40% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Paraneoplastic antigen Ma2 is a protein that in humans is encoded by the PNMA2 gene.
---------------------	--

GENE INFORMATION

Entrez Gene ID [10687](#)

Protein Refseq [NP_009188](#)

UniProt ID [Q5U5Z3](#)
