



# Anti-PNRC1 polyclonal antibody (CABT-BL2941)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                                     |
|--------------------|-----------------------------------------------------------------------------------------------------|
| Immunogen          | Recombinant protein corresponding to amino acids of human PNRC1.                                    |
| Isotype            | IgG                                                                                                 |
| Source/Host        | Rabbit                                                                                              |
| Species Reactivity | Human                                                                                               |
| Purification       | Antigen affinity purification                                                                       |
| Conjugate          | Unconjugated                                                                                        |
| Applications       | WB, IHC                                                                                             |
| Format             | Liquid                                                                                              |
| Buffer             | In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)                                                   |
| Preservative       | 0.02% Sodium Azide                                                                                  |
| Storage            | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

## BACKGROUND

|              |                                                                                                       |
|--------------|-------------------------------------------------------------------------------------------------------|
| Introduction | Proline-rich nuclear receptor coactivator 1 is a protein that in humans is encoded by the PNRC1 gene. |
|--------------|-------------------------------------------------------------------------------------------------------|

## GENE INFORMATION

Entrez Gene ID [10957](#)

---

Protein Refseq [NP\\_006804](#)

---

UniProt ID [Q12796](#)

---