



Anti-PODXL2 (aa 572-602) polyclonal antibody (CABT-BL2944)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 572-602 of Human Endoglycan.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Single-pass type I membrane protein
Format	Liquid
Buffer	99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	This gene is a member of the CD34 family of cell surface transmembrane proteins, which are characterized by an N-terminal extracellular mucin domain, globular and stalk domains, a single pass transmembrane region, and a charged cytoplasmic tail. The encoded protein is a ligand
---------------------	---

GENE INFORMATION

Entrez Gene ID	50512
----------------	-----------------------

Protein Refseq	NP_056535
----------------	---------------------------

UniProt ID	Q9NZ53
------------	------------------------
