



Anti-POTED (C-terminal) polyclonal antibody (CABT-BL2955)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------------------|---|
| Immunogen | KLH conjugated synthetic peptide selected from the C terminal region of human POTE. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human |
| Purification | Protein A purified |
| Conjugate | Unconjugated |
| Applications | IHC-P, ELISA |
| Cellular Localization | Cell membrane; Peripheral membrane protein. |
| Format | Liquid |
| Buffer | PBS |
| Preservative | 0.09% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

GENE INFORMATION

| | |
|-----------------------|---------------------------|
| Entrez Gene ID | 317754 |
| Protein Refseq | NP_778146 |

UniProt ID

[Q86YR6](#)
