



Anti-POTEG (internal region) polyclonal antibody (CABT-BL2957)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	The antiserum was produced against synthesized peptide derived from internal of human A26C2/3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	IHC, ELISA
Buffer	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C/1 year

GENE INFORMATION

Entrez Gene ID	404785
Protein Refseq	NP_001005356
UniProt ID	Q6S5H5