



Anti-PPFIA1 (aa 310-360) polyclonal Antibody (CABT-BL2958)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment: E AMAQKEDM EE RITTLEK RYL AAQREA TSVH DLNDK LENEI ANKD SMHRQT, corresponding to amino acids 310/360 of Human PPFIA1
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Mouse, Rat, Human, Fruit fly
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC
Cellular Localization	Cytoplasmic
Format	Liquid
Buffer	PBS
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	May regulate the disassembly of focal adhesions. May localize receptor-like tyrosine phosphatases type 2A at specific sites on the plasma membrane, possibly regulating their
---------------------	---

interaction with the extracellular environment and their association with substrates.

GENE INFORMATION

Entrez Gene ID	8500
----------------	----------------------

Protein Refseq	NP_003617
----------------	---------------------------

UniProt ID	B3KVS8
------------	------------------------
