



# Anti-PPFIA1 (aa 574-604) polyclonal antibody (CABT-BL2959)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein of mouse liprin-Alpha4(aa574- 604).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB, ICC
Preservative	None
Storage	-20°C, Avoid freeze / thaw cycles

## BACKGROUND

Introduction	May regulate the disassembly of focal adhesions. May localize receptor-like tyrosine phosphatases type 2A at specific sites on the plasma membrane, possibly regulating their interaction with the extracellular environment and their association with substrates.
--------------	---

## GENE INFORMATION

Entrez Gene ID	<a href="#">8500</a>
Protein Refseq	<a href="#">NP_003617</a>

UniProt ID

[B3KVS8](#)

---