



# Anti-PPFIBP1 (aa 37-86) polyclonal antibody (CABT-BL2960)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: SPTSPFMGSL RALHLVEDLR GLEMMETDE KEGLRCQIPD STAETLVEWL, corresponding to amino acids 37-86 of Human PPFIBP1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

Introduction	May regulate the disassembly of focal adhesions. Did not bind receptor-like tyrosine phosphatases type 2A.
--------------	--

## GENE INFORMATION

Entrez Gene ID [8496](#)

---

Protein Refseq [NP\\_001185844](#)

---

UniProt ID [A0A024RB02](#)

---