



Anti-PPIAL3 (full length) polyclonal antibody (CABT-BL2962)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full length protein (Human).
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	IgY fraction
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [375113](#)