



Anti-PPP2R5A (aa 473-486) polyclonal antibody (CABT-BL2976)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: AYNMHSILSNTSAE, corresponding to C terminal amino acids 473-486 of Human PPP2R5A.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm. Nucleus. Chromosome; centromere. From mitotic prophase to metaphase, localizes at the inner centromere between a pair of sister kinetochores. Decreased expression at the onset of anaphase.
Format	Liquid
Buffer	0.5% BSA, Tris buffered saline, pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

GENE INFORMATION

Entrez Gene ID	5525
----------------	----------------------

Protein Refseq	NP_001186685
----------------	------------------------------

UniProt ID	Q15172
------------	------------------------
