



Anti-PPP2R5A (aa 25-75) polyclonal antibody (CABT-BL2977)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from a region between residue 25 and 75 of human PPP2R5A, NP_006234.1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Cytoplasm. Nucleus. Chromosome; centromere. From mitotic prophase to metaphase, localizes at the inner centromere between a pair of sister kinetochores. Decreased expression at the onset of anaphase.
Format	Liquid
Size	100 µl
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

GENE INFORMATION

Entrez Gene ID	5525
----------------	----------------------

Protein Refseq	NP_001186685
----------------	------------------------------

UniProt ID	Q15172
------------	------------------------
