



Anti-PPP2R5A polyclonal antibody (CABT-BL2979)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 108-157 (VSTNRGVIVE SAYSDIVKMI SANIFRTLPP SDNPFDPEE DEPTLEASWP) of human PPP2R5A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasm. Nucleus. Chromosome; centromere. From mitotic prophase to metaphase, localizes at the inner centromere between a pair of sister kinetochores. Decreased expression at the onset of anaphase.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

GENE INFORMATION

Entrez Gene ID	5525
----------------	----------------------

Protein Refseq	NP_001186685
----------------	------------------------------

UniProt ID	Q15172
------------	------------------------
