



Anti-PPP2R5C (aa 481-490) polyclonal antibody (CABT-BL2982)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein: QAQKDPKKDR, corresponding to amino acids 481-490 of Human PPP2R5C
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus. Chromosome; centromere.
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.
---------------------	---

GENE INFORMATION

Entrez Gene ID [5527](#)

Protein Refseq [NP_001155197](#)

UniProt ID [Q13362](#)
