



Anti-PPP4R2 polyclonal antibody (CABT-BL2987)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residue 370 and the C-terminus (residue 417) of Human PPP4R2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Cellular Localization	Centrosome.
Format	Liquid
Size	100 µl
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction PPP4R2 (protein phosphatase 4, regulatory subunit 2) is a protein-coding gene. Diseases

associated with PPP4R2 include hepatitis, and neuronitis. GO annotations related to this gene include protein phosphatase type 4 regulator activity and protein binding, bridging.

GENE INFORMATION

Entrez Gene ID	151987
----------------	------------------------

Protein Refseq	NP_777567
----------------	---------------------------

UniProt ID	Q9NY27
------------	------------------------
