



Anti-PQLC2 (aa 72-121) polyclonal antibody (CABT-BL2991)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from within residues: LFLMGMACAT PLLSAAGPVA APREAFRGRA LLSVESGSKP FTRQEVIGFV, corresponding to internal sequence amino acids 72-121 of Human PQLC2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Integral to membrane
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The specific function of this protein remains unknown. There are 3 named isoforms.

GENE INFORMATION

Entrez Gene ID [54896](#)

Protein Refseq [NP_001035214](#)

UniProt ID [A0A024RA97](#)
