



Anti-PRAMEF3 (aa 427-457) polyclonal antibody (CABT-BL2993)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 427-457 of Human PRAMEF3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	DISCONTINUED: This record has been withdrawn by NCBI. This sequence is no longer present in the current reference genome assembly.
--------------	--

GENE INFORMATION

Entrez Gene ID [401940](#)

Protein Refseq [NP_001013714](#)
