



Anti-PRD (aa 1-50) polyclonal antibody (CABT-BL2994)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: MTVTAFAAAM HRPFFNGYST MQDMNSGQGR VNQLGGVFIN GRPLPNNIRL, corresponding to N terminal amino acids 1-50 of Fruit fly (Drosophila melanogaster) Segmentation protein paired.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, 1X PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [34629](#)

Protein Refseq [NP_723721](#)

UniProt ID [P06601](#)
