



Anti-PRDM11 (aa 499-511) polyclonal antibody (CABT-BL2995)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: LSEGRVRSGLCGG, corresponding to C terminal amino acids 499-511 of Human PRDM11
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	0.5% BSA, 0.5mg/ml Tris buffered saline. pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	PRDM11 (PR domain containing 11) is a protein-coding gene. An important paralog of this gene is PRDM9.
---------------------	--

GENE INFORMATION

Entrez Gene ID [56981](#)

Protein Refseq [NP_064614](#)

UniProt ID [A0A087WWZ6](#)
