



Anti-PRDM12 (aa 161-210) polyclonal antibody (CABT-BL2996)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 161-210 (MTYIKCARNE QEQNLEVVQI GTSIFYKAIE MIPPDQELLV WYGNSHNTFL) of Mouse PRDM12.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, 98% PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	PRDM12 (PR domain containing 12) may be involved in transcriptional regulation.
---------------------	---

GENE INFORMATION

Entrez Gene ID [381359](#)

Protein Refseq [NP_001116834](#)

UniProt ID [A2AJ77](#)
