



# Anti-PRDM13 (aa 649-698) polyclonal antibody (CABT-BL2999)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                           |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | A region within synthetic peptide: VKSRHPGQSL LAKAGDGPGA EPGYPPEPGD PKSDSDVDV CFTDDQSDPE, corresponding to C terminal amino acids 649-698 of Human PRDM13 |
| <b>Isotype</b>               | IgG                                                                                                                                                       |
| <b>Source/Host</b>           | Rabbit                                                                                                                                                    |
| <b>Species Reactivity</b>    | Human                                                                                                                                                     |
| <b>Purification</b>          | Protein A purified                                                                                                                                        |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                              |
| <b>Applications</b>          | WB, ELISA                                                                                                                                                 |
| <b>Cellular Localization</b> | Nuclear                                                                                                                                                   |
| <b>Format</b>                | Liquid                                                                                                                                                    |
| <b>Buffer</b>                | 2% Sucrose, PBS                                                                                                                                           |
| <b>Preservative</b>          | None                                                                                                                                                      |
| <b>Storage</b>               | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.                                                                     |

## GENE INFORMATION

Entrez Gene ID [59336](#)

Protein Refseq [NP\\_067633](#)

---

UniProt ID [Q9H4Q3](#)

---