



Anti-PRE3 (aa 61-160) polyclonal antibody (CABT-BL3000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Vector coding for a partial recombinant fusion protein, corresponding to amino acids 61-160 of <i>Saccharomyces cerevisiae</i> Pre3. Target sequence used to make the antibody: WCCRSGSAAD TQAIADIVQY HLELYTSQYG TPSTETAASV FKELCYENKD NLTAGIIVAG YDDKNKGEVY TIPLGGSVHK
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Entrez Gene ID [853456](#)

Protein Refseq [NP_012533](#)

UniProt ID [P38624](#)
