



Anti-PRE7 (aa 11-110) polyclonal antibody (CABT-BL3001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Vector coding for a partial recombinant fusion protein, corresponding to amino acids 11-110 of <i>Saccharomyces cerevisiae</i> PRE7. Target sequence used to make the antibody:- YSFDPVGSYEREQCRAAGAAASLIMPFLDNQVNFKNQY EPGTNGKVKKPLKYLSVEEVIKLVRDSFTSATERHIQVGD GLEILI
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [852239](#)

Protein Refseq [NP_009512](#)

UniProt ID [P23724](#)
