



Anti-prg-2 polyclonal antibody (CABT-BL3005)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 700 to the C-terminus of <i>Caenorhabditis elegans</i> PRG2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>C. elegans</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Size	100 µg
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [182072](#)

Protein Refseq [NP_500994](#)

UniProt ID [Q17567](#)
