



Anti-PRLH polyclonal antibody (CABT-BL3020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Buffer	PBS with 0.09% (W/V) sodium azide
Preservative	0.09% Sodium Azide
Storage	PRLH Antibody (Center) can be refrigerated at 2-8°C for up to 6 months. For long term storage, place the at -2°C.

GENE INFORMATION

Entrez Gene ID	623503
Protein Refseq	NP_001095117
UniProt ID	G3UWC3