



Anti-PROS26.4 (aa 8-107) polyclonal antibody (CABT-BL3025)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: QGGAGEKKDDKDKKKKYEPPIPTRVGKKKRRAKGPDAAMK LPQVTPHTRCRLKLLKLERIKDYLMMEDEFIRNQERLKPQ DEKNEEERSKVDDLGRGTPMS, corresponding to amino acids 8/107 of Drosophila melanogaster RPT2.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Fruit fly
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic and Nuclear
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [42828](#)

Protein Refseq [NP_524469](#)

UniProt ID [P48601](#)
