



Anti-PRR13 (aa 35-84) polyclonal antibody (CABT-BL3027)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 35-84 (PFPPGPCPPP PGAPHGNAF PPGGPPHPVP QPGYPGCQPL GPYPPPYPPP) of Human PRR13.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	54458
Protein Refseq	NP_001005354

UniProt ID

[Q9NZ81](#)
