



Anti-PRR15 (aa 79-128) polyclonal antibody (CABT-BL3028)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 79-128 (GDKSGSSRRN LKISRSGRFK EKRKVRATLL PEAGRSPEEA GFPDPHEDK) of Human PRR15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	222171
Protein Refseq	NP_787083

UniProt ID

[A4D1A1](#)
