



Anti-PRRC2C (aa 2250-2300) polyclonal antibody (CABT-BL3031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within amino acids 2250-2300 of Human PRRC2C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	99% Tris citrate/phosphateNote: pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	23215
Protein Refseq	NP_055987

UniProt ID

[Q9Y520](#)
