



Anti-PRRG3 (aa 38-87) polyclonal antibody (CABT-BL3032)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 38-87 (EEICSYEEVK EVFENKEKTM EFWKGYPNAV YSVRDPSQSS DAMYVVVPLL) of Human PRRG3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes a protein which contains a vitamin K-dependent carboxylation/gamma-carboxyglutamic domain. The encoded protein is a member of a family of vitamin K-dependent
---------------------	---

transmembrane proteins which contain a glutamate-rich extracellular domain.

GENE INFORMATION

Entrez Gene ID	79057
----------------	-----------------------

Protein Refseq	NP_076987
----------------	---------------------------

UniProt ID	Q9BZD7
------------	------------------------
