



Anti-PRUNE2 polyclonal antibody (CABT-BL3037)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human PRUNE2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	PRUNE2 (prune homolog 2 (Drosophila)) is a protein-coding gene. Diseases associated with PRUNE2 include leiomyosarcoma, and gastrointestinal stromal tumor. GO annotations related to this gene include pyrophosphatase activity and metal ion binding. An important paralog of this gene is BNIP2.
---------------------	---

GENE INFORMATION

Entrez Gene ID [158471](#)

Protein Refseq [NP_056040](#)

UniProt ID [Q8WUY3](#)
