



Anti-psbA (aa 1-17) polyclonal antibody (CABT-BL3041)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: (M)TAILERRESESLWGRFC, corresponding to N terminal amino acids 1-17 of Arabidopsis thaliana psbA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana, Tobacco
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB Recommended dilution: WB: 1:2500---1:5000 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Cellular Localization	Plastid, chloroplast thylakoid membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	Preservative: 0.01% Sodium Azide
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [844802](#)

Protein Refseq [NP_051039](#)

UniProt ID [P83755](#)
