



Anti-psbD (aa 230-240) polyclonal antibody (CABT-BL3042)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: ANTFRANPTQ, corresponding to amino acids 230-240 of Arabidopsis thaliana psbD
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana, Rice
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Plastid, chloroplast thylakoid membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	Whole serum
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [844775](#)

Protein Refseq [NP_051054](#)

UniProt ID [P56761](#)
