



# Anti-PSMG3 polyclonal antibody (CABT-BL3044)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human PSMG3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Introduction	PSMG3 is a chaperone protein which promotes assembly of the 20S proteasome. This protein may cooperate with PSMG1-PSMG2 heterodimers to orchestrate the correct assembly of proteasomes.
--------------	--

# GENE INFORMATION

Entrez Gene ID [84262](#)

---

Protein Refseq [NP\\_001127812](#)

---

UniProt ID [A0A024R806](#)

---