



Anti-PSORS1C2 (internal region) polyclonal antibody (CABT-BL3045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	The antiserum was produced against synthesized peptide derived from internal of human SPR1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C/1 year

GENE INFORMATION

Entrez Gene ID	170680
Protein Refseq	NP_054788
UniProt ID	Q9UIG4