



Anti-PSTK (aa 144-193) polyclonal antibody (CABT-BL3046)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 144 -193 (RPLFLVLDDNFYYQSMRYEVYQLARKYSLGFCQLFLDCP LETCLQRNGQR) of Human PSTK
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [118672](#)

Protein Refseq [NP_699167](#)

UniProt ID

[Q8IV42](#)
