



Anti-PTGIS (full length) polyclonal antibody (CABT-BL3049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length human PTGIS protein
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Endoplasmic reticulum membrane.
Format	Liquid
Buffer	PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Catalyzes the isomerization of prostaglandin H2 to prostacyclin (= prostaglandin I2).
--------------	---

GENE INFORMATION

Entrez Gene ID [5740](#)

Protein Refseq [NP_000952](#)

UniProt ID [Q16647](#)
