



Anti-PTGIS (aa 475-490) polyclonal antibody (CABT-BL3050)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: PEFDLSTRYGFLMQPE, corresponding to amino acids 475-490 of Mouse Prostaglandin I synthase (PTGIS).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Sheep, Cow, Human
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, IHC-Fr, WB
Cellular Localization	Endoplasmic reticulum membrane.
Format	Liquid
Size	100 µl
Buffer	Preservative: None
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Catalyzes the isomerization of prostaglandin H2 to prostacyclin (= prostaglandin I2).
---------------------	---

GENE INFORMATION

Entrez Gene ID [282021](#)

Protein Refseq [NP_776869](#)

UniProt ID [Q29626](#)
