



Anti-PTH2 (aa 62-100) polyclonal antibody (CABT-BL3053)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	synthetic human TIP 39 (aa 62-100) bTG conjugated (SLALADDAAFRERARLLAALERRHWLNSYMHKLLVLDAP)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA, RIA
Format	Serum
Preservative	None
Storage	2°C-8°C (lyophilized); -20°C (dissolved) Repeated thawing and freezing must be avoided

BACKGROUND

Introduction	TIP39 is related to parathyroid hormone (PTH; MIM 168450) and PTH-related protein (PTHrP; MIM 168470) and is a ligand;for PTH receptor-2 (PTHR2; MIM 601469) (John et al., 2002 (PubMed 11861531)).
---------------------	---

GENE INFORMATION

Entrez Gene ID	113091
-----------------------	------------------------

Protein Refseq [NP_848544](#)

UniProt ID [Q96A98](#)
