



Anti-PTP2 (aa 280-379) polyclonal antibody (CABT-BL3055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: ILPNWLKFCSVKENEKVILKKLFNNFETLENFEMQRLEKC LKFKKKPLHQKQLSQKQRGPQSTDDSKLYSLTSLQRQYKS SLKSNIQKNQKCLKLIIPKNN, corresponding to amino acids 280/379 of <i>S. cerevisiae</i> PTP2
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cytoplasmic
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [854383](#)

Protein Refseq [NP_014851](#)

UniProt ID [P29461](#)
